

# Two Dimensional Homonuclear Correlation Spectroscopy

## Gradient COSY

William D. Wheeler, Ph.D.  
Department of Chemistry  
University of Wyoming

April 16, 1999  
Revised September 22, 1999

## INTRODUCTION

### Correlation Spectroscopy


Correlation experiments reveal connections between spin coupled nuclei. Spin, or scalar, coupling between nuclei is established by the appearance of cross peaks (peaks off the diagonal) in the two-dimensional spectrum. Since spin coupled nuclei are usually separated by two or three bonds, the correlation experiment, by revealing the connectivity within the compound, is often enough to establish the chemical structure. Resonances which are coupled through four and five bonds will occasionally show cross peaks also. Although most correlation experiments are designed to emphasize short range coupling, there are experiments designed to emphasize long range correlations as well.

Two-dimensional spectra can be acquired in either of two modes. In the magnitude, or absolute value mode, only the "real" part of the data is collected since the spectrum is not phased. In the phase sensitive mode, both the "real" and "imaginary" parts of the data are acquired so that the spectrum can be phased. An advantage of a phase sensitive spectrum is that the peaks are much narrower and the resolution greater. A disadvantage of a phase sensitive spectrum is that twice as much data must be recorded in order to extract the phase information. Thus, it takes twice the time to record a spectrum in phase sensitive mode.


### Gradient COSY

The experiment described here is the simplest of two-dimensional magnetic field gradient NMR experiments. Magnetic field gradients are sometimes called pulsed field gradients, pfg. Normal COSY spectra require a minimum of four scans to distinguish the sign of a frequency in the second ( $F_1$ ) dimension and select the desired coherence. If you require the peaks to have the correct shape (phase sensitive mode), you must collect eight scans. By introducing magnetic field gradients into the experiment, this can all be accomplished in one scan. Gradient methods can significantly shorten the time it takes to acquire a two-dimensional spectrum, provided you have a concentrated enough sample. An additional benefit of gradients is a significantly reduced  $t_1$  noise. Gradient COSY spectra can be acquired in both magnitude and phase sensitive mode.

The pulse sequence for gradient-COSY is shown below.




The gradient COSY spectrum of Menthol is shown on a following page. The spectrum required about 11 minutes to collect.



### References.

- 1) **150** and More Basic NMR Experiments, A practical Course, S. Braun, H.-O. Kalinowski, S. Berger, Wiley-VCH, 1998, pages 477-480.



```

Current Data Parameters
NAME      Menthol_COSY
EXPNO     5
PROCNO    1

F2 - Acquisition Parameter
Date_     990223
Time      13.09
INSTRUM   spect
PROBHD    5 mm QNP 1H
PULPROG   cosygp
TD         1024
SOLVENT   CDCl3
NS         1
DS         16
SWH        2003.205 Hz
FIDRES     1.956255 Hz
AQ         0.2556404 se
RG         8
DW         249.600 us
DE         6.00 us
TE         300.0 K
d0         0.00000300 se
d13        0.00000300 se
D1         1.00000000 se
P1         7.40 us
SF01       400.1308094 MHz
NUC1       1H
P11        -3.00 dB
P16        1000.00 us
GPX1       0.00 %
GPY1       0.00 %
GPZ1       10.00 %
GNAM1      sine 100
D16        0.00025000 se
P0         7.40 us
GPX2       0.00 %
GPY2       0.00 %
GPZ2       10.00 %
GNAM2      sine 100
INO        0.00049920 se

F1 - Acquisition parameter
ND0        1
TD         256
SF01       400.1308 MHz
FIDRES     7.825020 Hz
SW         5.006 ppm

F2 - Processing parameters
SI         512
SF         400.1300092 MHz
WDW        SINE
SSB        0
LB         0.00 Hz
GB         0
PC         1.00

F1 - Processing parameters
SI         512
MC2        GF
SF         400.1300092 MHz
WDW        SINE
SSB        0
LB         0.00 Hz
GB         0

2D NMR plot parameters
CX2        15.00 cm
CX1        15.00 cm
F2PLO     3.604 ppm
F2LO      1441.94 Hz
F2PHI     0.396 ppm
F2HI      158.63 Hz
F1PLO     3.604 ppm
F1LO      1441.94 Hz
F1PHI     0.396 ppm
F1HI      158.63 Hz
F2PPMCM   0.21381 ppm
F2HZCM    85.55355 Hz
F1PPMCM   0.21381 ppm
F1HZCM    85.55355 Hz
    
```

Gradient COSY of Menthol

## EXPERIMENTAL SETUP

The concentration of the sample should be high enough that you can record a high signal to noise  $^1\text{H}$  NMR spectrum in 16 - 32 scans.

### Record a $^1\text{H}$ spectrum.

NEW (or EDC)	Create a new data set for your sample.
NAME name	Data set name.
EXPNO 1	Experiment number (must be 1).
PROCNO 1	Process data number (must be 1).
[SAVE]	Save the data set.
RPAR +proton all	Read in the standard proton parameters.
NS, etc.	Adjust NS and other parameters as necessary.
RGA, ZG	Acquire some data.
FT, APK, ref	Fourier transform, phase and reference the spectrum.
Zoom	Zoom in on the region of interest. The expanded region need not contain a reference peak.
[sw-sfo1]	Set the sweep width and spectrometer frequency to cover the zoomed region. Reduce TD if the acquisition time (AQ) is unnecessarily large.
RGA, ZG	Acquire a spectrum of the zoomed region.
FT, etc.	Fourier transform, phase etc.

### Setting the conditions for $^1\text{H}$ - $^1\text{H}$ correlation.

The following AU program sets up all of the parameters for gradient cosy.

XAU su_gcosy	Set up gradient cosy (g cosy).
[OK]	Answer "Delete `meta.ext` files ?" with [OK].

The AU program sets the following file parameters. Use EDC to display the results.

File parameter		
EXPNO	1	Experiment number 1, for 1-D proton.
	5	Experiment number 5, for g cosy.

The AU program sets the following acquisition parameters. Use EDA to display the results.

Acquisition Parameters		
<b>Time domain 2 (F2)</b>		
PULPROG	cosygp	PULse PROGram for m cosy.
TD	1024	Time Domain points.
NS	1	Number of Scans (should be 1, 2, 4, 8 etc.).
DS	16	number of Dummy Scans for 1st row.
D1	2 sec	relaxation Delay.
GPZ1	10.0	strength of the first Gradient Pulse along Z.
GPZ2	10.0	strength of the second Gradient Pulse along Z.
<b>Time domain 1 (F1)</b>		
TD	256	Time Domain points.
ND0	1	Number of D0 periods per cycle.
IN0	1 / SW	Increment in t1 (calculated).
SW		SW of <sup>1</sup> H spectrum (same as for F2).
SFO1		<sup>1</sup> H frequency (same as for F2).

ASED

Check that acquisition parameters are set correctly. A brief description of other parameters, not described above, is given in the pulse program at the end of this document.

[SAVE]

Save the acquisition parameters.

EXPT

Calculates the approximate length of time to do the experiment. This may help you to decide if you want to collect more or fewer slices or scans or points etc. The Menthol spectrum required 11 minutes.

## DATA ACQUISITION

[Spin on/off]

TURN THE SPINNER OFF ( press the Spin On/Off button on the BOSS keyboard).

ZG

Start the acquisition.

## DATA PROCESSING

The AU program sets the following processing parameters. Use EDP to check or modify the values.

Processing Parameters		
<b>Time/frequency domain 2 (F2)</b>		
SI	512	the SIze in F2 (zero fill rows).
WDW	SINE	Sine multiplication.
SSB	0	Unshifted sine bell.
PH_mod	NO	No phase correction.
PKNL	TRUE	Required with digital filter.
BC_mod	QUAD	Background correct quadrature data.
<b>Time/frequency domain 1 (F1)</b>		
SI	512	the SIze in F1 (zero fill columns).
WDW	SINE	Sine multiplication.
SSB	0	Unshifted sine bell.
PH_mod	MC	Magnitude calculation.
BC_mod	NO	No background correction.
MC2	QF	Forward quadrature FFT.
OFFSET		Frequency offset (same as for F2).
SF		Spectrometer frequency (same as for F2).

XFB Background correct, window, zero fill, Fourier transform, phase and reference. The whole kaboodle in one command. It is OK to execute this command on a partial data set, during acquisition.

## 2-D CONTOUR DISPLAY

### X-Y Expansion the spectrum.

[Limits] Set the limits of the plot region. NOTE: For homo-nuclear correlation experiments, (COSY, NOESY, etc.) set F1LO = F2LO and F1HI = F2HI so the display will be symmetric and the diagonal peaks will appear on the diagonal line.

[PlotReg] Forces XWinNMR to display the plot region.

### Setting the intensity scale.

[DefPlot] Sets the intensity scale for plotting to be the same as the currently displayed intensity scale.

[contours] Displays the equi-intensity contours of the data.

[intensities] Displays ranges of intensity as a color map.



<b>Additional Processing Parameters</b>		
<b>Frequency domain</b>		
F1LO	bbb	Set to the value for the bottom of the spectrum (ppm).
F1HI	ttt	Set to the value for the top of the spectrum (ppm).
F2LO	lll	Set to the value for the left side of the spectrum (ppm).
F2HI	rrr	Set to the value for the right side of the spectrum (ppm).

<b>Additional Processing Commands</b>	
<b>Frequency domain</b>	
ABS1	Automatic Baseline Subtraction for F1 (uses F1LO, F1HI).
ABS2	Automatic Baseline Subtraction for F2 (uses F2LO, F2HI).
SYM	Symmetrize the spectrum. This “averages” the upper and lower triangles of the matrix. CAUTION! This is an easy way to introduce artifacts into your spectrum.

## Pulse Program for Gradient COSY

```

;cosygp
;avance-version
;2D homonuclear shift correlation
;using gradient pulses for selection

#include <Avance.incl>
#include <Grad.incl>

"d0=3u"
"d13=3u"

1 ze
2 d1
3 p1 ph1
  d0
  50u UNBLKGRAD
  p16:gp1
  d16
  p0 ph2
  d13
  p16:gp2
  d16
  4u BLKGRAD
  go=2 ph31
  d1 wr #0 if #0 id0 zd
  lo to 3 times tdl
exit

ph1=0 2
ph2=0 0 2 2
ph31=0 2

;p11 : f1 channel - power level for pulse (default)
;p0  : f1 channel - 20 to 90 degree high power pulse
;p1  : f1 channel - 90 degree high power pulse
;p16: homospoil/gradient pulse
;d0  : incremented delay (2D) [3 usec]
;d1  : relaxation delay; 1-5 * T1
;d13: short delay [3 usec]
;d16: delay for homospoil/gradient recovery
;in0: 1/(1 * SW) = 2 * DW
;nd0: 1
;NS: 1 * n
;DS: 16
;td1: number of experiments
;MC2: QF
;use gradient ratio:      gp 1 : gp 2
;                          10 : 10

```